

Hammertoe

PEEK Fixation System



Guided PEEK solution for hammertoe, claw toe, and mallet toe procedures




Hammertoe

PEEK Fixation System

- **PEEK implant**
offers radiolucency and mechanical properties comparable to bone^{1,2}
- **Cannulated system**
facilitates streamlined insertion and allows for accurate implant placement
- **Sterile-packed**
implant and instrument kit for convenient delivery and clinical efficiency



Ordering Information

Small		Medium		Large	
					
Ref#	Description	Size	Drill	Tap	Guidewire
SK31	Hammertoe PEEK Fixation System	Small	2.75mm	2.5mm	1.1mm Double Trocar
SK32	Hammertoe PEEK Fixation System	Medium	3.15mm	3.1mm	1.4mm Double Trocar
SK33	Hammertoe PEEK Fixation System	Large	3.45mm	3.7mm	1.6mm Double Trocar

Before use of the system, the surgeon should refer to the appropriate instructions for use and surgical technique for complete warnings, precautions, indications, contraindications, and adverse events. Risks include, but are not limited to: infection, pain, discomfort from the presence of the implant, loosening of the implant, and loss of correction with nonunion or malunion. If any of these occur, additional treatments may be needed. Additional information about risks, warnings, and instructions is available at [Lapiplasty.com/surgeons/labeling](https://lapiplasty.com/surgeons/labeling).

To learn more, visit
Lapiplasty.com



TREACE
Medical Concepts, Inc.